


AMICO 3-GAS ARTWALL SYSTEM

(M/N: W-ARTWALL-03R)

DRAWING # 45




MOUNTING PLATE
CAN ACCOMMODATE
35"x25" FRAMED PRINTS
(Provided by customer, or
Amico at additional charge).

TYPE:
QUANTITY:

SYSTEM DETAILS	
SYMBOL	DESCRIPTION
O	GAS, SWIVEL TYPE, OXYGEN
A	GAS, SWIVEL TYPE, MED AIR
V	GAS, SWIVEL TYPE, VACUUM
SL	VACUUM SLIDE, FLIP DOWN

- FEATURES:
- RECESSED ENCLOSURE (Ga. #16 PAINTED STEEL)
 - SWIVEL GAS OUTLETS
 - FLOWMETERS & REGULATORS CAN BE ATTACHED TO SWIVEL OUTLETS
 - UPWARD PICTURE FRAME (DIAGONAL LEFT OR RIGHT)
 - ACCEPTS FRAMED WORK
 - GASES MANIFOLD ON SITE BY DIV. 25 CONTRATOR



IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE _____

DATE _____

PHONE NO. _____

Amico Corporation
 85 Fulton Hwy, Ontario
 Richmond Hill, Canada
 L4B 2N4, CANADA
 Toll-Free: 1-877-462-6428 (T)
 Tel: (905) 764-0980
 Fax: (905) 764-0982
 www.amico.com

HOSPITAL LOCATION: _____
 HOSPITAL LOCATION: _____
 QTY: 1 () TYPE _____ UNITS AS SHOWN / () TYPE _____ UNITS OPPOSITE

A. NURSE CALL MFGR.: _____ MODEL NO.: _____
 B. MEDICAL GAS MFGR.: _____ TYPE CONNECTION: _____
 C. FINISH: _____ CEILING HEIGHT: _____
 DRWG. NO. PRICEBOOK-49
 DRAWN BY: _____
 CHECKED BY: _____
 REV. NO.: _____
 DATE: _____