

# REGAL SERIES RECESSED 19" FLATWALL EXTENDED LENGTH GENERAL CARE - W/STUD SINGLE SIDE


(M/N: FW1907-RC-GEN-P)

DRAWING # 08

TYPE:  
QUANTITY:


SYSTEM DETAILS		
SYMBOL	QTY.	DESCRIPTION
FW	1	FLATWALL, 19" WIDTH
O	3	GAS, OXYGEN
A	2	GAS, MED AIR
V	2	GAS, VACUUM
E2	5	RECEPTACLE - DUPLEX RED
N2	4	RECEPTACLE - DUPLEX IVORY
BR	1	RECEPTACLE - BED ONLY
SN1	1	SWITCH, SPST
NC	1	PROVISION - NURSE CALL 3G
CB	1	PROVISION - CODE BLUE 1G
DA	3	PROVISION - DATA 1G
TL	1	PROVISION - TELEPHONE 1G
SC	1	37 PIN BED CONNECTER
MC	1	MONITOR CHANNEL, MOUNTED
	1	STUD ASSEMBLY

NOTE:  
-STUD ASSEMBLY PROVIDED BY AMICO



**IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.**

APPROVAL SIGNATURE \_\_\_\_\_ DATE \_\_\_\_\_ PHONE NO. \_\_\_\_\_



**AMICO Corporation**  
85 Fulton Way  
Richmond Hill, Ontario  
L4B 2M4, CANADA  
Tel: (905) 1-866-440-4886  
Fax: (905) 764-0892  
www.amico.com

HOSPITAL LOCATION: \_\_\_\_\_  
HOSPITAL LOCATION: \_\_\_\_\_  
QTY: (A) TYPE \_\_\_\_\_ UNITS AS SHOWN / (B) TYPE \_\_\_\_\_ UNITS OPPOSITE

A. NURSE CALL MFR: \_\_\_\_\_ MODEL NO.: \_\_\_\_\_  
B. MEDICAL GAS MFR: \_\_\_\_\_ TYPE CONNECTION: \_\_\_\_\_  
C. FINISH: \_\_\_\_\_ CEILING HEIGHT: \_\_\_\_\_  
DRWG. NO. 01  
PRICEBOOK-07  
DRAWN BY:  
CHECKED BY:  
REV. NO.: 00  
DATE: