


REGAL SERIES SURFACE MOUNTED 19" FLATWALL


DRAWING # 14

FIVE EIGHTHS LENGTH GENERAL CARE
(M/N: FW1906-SM-GEN-P)



TYPE: _____
QUANTITY: _____

SYSTEM DETAILS		
SYMBOL	QTY.	DESCRIPTION
FW	1	FLATWALL, 19" WIDTH
O	1	GAS, OXYGEN
A	1	GAS, MED AIR
V	1	GAS, VACUUM
N2	2	RECEPTACLE - DUPLEX IVORY
NC	1	PROVISION - NURSE CALL 3G
SLM	1	VACUUM SLIDE, TRACK MTD.
VVH	1	I.V.HOOK, TRACK MTD.



IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE _____

DATE _____ PHONE NO. _____

AMICO Corporation
85 Fulton Hwy,
Rehoboth Hill, Ontario
L4B 2N4, CANADA
Tel: (905) 877-462-6426 (T)
Toll-Free: 1-877-462-6426 (T)
Fax: (905) 764-0800
www.amico.com

HOSPITAL	HOSPITAL
LOCATION	LOCATION
QTY.	1
(A) TYPE	UNITS AS SHOWN / (B) TYPE
	UNITS OPPOSITE

A. NURSE CALL MFGR.:	MODEL NO.:	DRWG. NO.:
B. MEDICAL GAS MFGR.:	TYPE CONNECTION:	PRICEBOOK-17
C. FINISH:	CEILING HEIGHT:	DRAWN BY:
		CHECKED BY:
		REV. NO.:
		DATE: