


# REGAL SERIES SURFACE MOUNTED 38" FLATWALL


DRAWING # 24

FIVE EIGHTHS LENGTH GENERAL CARE  
(M/N: FW3806-SM-GEN-S)



TYPE:  
QUANTITY:

SYSTEM DETAILS	
SYMBOL	DESCRIPTION
FW	FLATWALL, 38" WIDTH
O	GAS, OXYGEN
V	GAS, VACUUM
N2	RECEPTACLE - DUPLEX IVORY
NC	PROVISION - NURSE CALL 3G



IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE \_\_\_\_\_ DATE \_\_\_\_\_ PHONE NO. \_\_\_\_\_

**AMICO Corporation**  
85 Fulton Hwy,  
Richmond Hill, Ontario  
L4B 2N4, CANADA  
Tel: (905) 764-0800  
Fax: (905) 764-0882  
www.amico.com

HOSPITAL	HOSPITAL	MODEL NO.:	DRWG. NO.
LOCATION	LOCATION	A. NURSE CALL MFGR.:	PRICEBOOK-27
QTY.:	TYPE	B. MEDICAL GAS MFGR.:	DRAWN BY:
1	UNITS AS SHOWN / ( ) TYPE UNITS OPPOSITE	C. FINISH:	CHECKED BY:
		CEILING HEIGHT:	REV. NO.:
			DATE: