



# REGAL SERIES 9" POWER COLUMN

FULL LENGTH  
(M/N: PC09FL-1CU)


DRAWING # 26



TYPE:  
QUANTITY:

SYSTEM DETAILS

SYMBOL	QTY.	DESCRIPTION
PC	1	POWER COLUMN, 9"x9"
FM	1	FLOOR MOUNTED SS BASE
VT	8	VERTICAL EQUIPMENT TRACK
O	2	GAS, OXYGEN
A	2	GAS, MED AIR
V	2	GAS, VACUUM
N2	4	RECEPTACLE - DUPLEX IVORY
E2	4	RECEPTACLE - DUPLEX RED
NC	2	PROVISION - NURSE CALL 3G
CB	1	PROVISION - CODE BLUE 1G
DA	2	PROVISION - DATA 1G



IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE \_\_\_\_\_

DATE \_\_\_\_\_

PHONE NO. \_\_\_\_\_

**AMICO Corporation**  
85 Fulton Hwy,  
Richmond Hill, Ontario  
L4B 2N4, CANADA  
Tel: (905) 764-0980  
Fax: (905) 764-0982  
www.amico.com

HOSPITAL	HOSPITAL	A. NURSE CALL MFR:	MODEL NO.:	DRWG. NO.
LOCATION	LOCATION	B. MEDICAL GAS MFR.:	TYPE CONNECTION:	PRICEBOOK-29
QTY. 1	( ) TYPE UNITS AS SHOWN / ( ) TYPE UNITS OPPOSITE	C. FINISH:	CEILING HEIGHT:	DRAWN BY:
				CHECKED BY:
				REV. NO.:
				DATE: