



REGAL SERIES 28" POWER COLUMN

DRAWING # 29

FULL LENGTH
(M/N: PC28FL-ICU)



IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.


APPROVAL - PRINT AND SIGNATURE _____

DATE _____ PHONE NO. _____

TYPE: _____
QUANTITY: _____

SYSTEM DETAILS	
SYMBOL	DESCRIPTION
PC	1 POWER COLUMN, 28"x9"
FM	1 FLOOR MOUNTED SS BASE
VT	8 VERTICAL EQUIPMENT TRACK
O	4 GAS, OXYGEN
A	1 GAS, MED AIR
V	2 GAS, VACUUM
E2	7 RECEPTACLE - DUPLEX RED
N2	6 RECEPTACLE - DUPLEX IVORY
NC	1 PROVISION - NURSE CALL 3G
CB	2 PROVISION - CODE BLUE 1G
DA	1 PROVISION - DATA 1G
DT	1 PROVISION - DATA QUAD 2G
MO	1 PROVISION - MONITOR OUT 2G
NL	2 NIGHT LIGHT
MC	1 MONITOR CHANNEL, MOUNTED
AM	1 ANEROID MANOMETER w/ BSKT

** 38 1/2" POWER COLUMN ALSO AVAILABLE**



AMICO Corporation
85 Fulton Hwy, Ontario
Rehoboth Hill, CANADA
L4B 2N4, CANADA
Tel: (905) 764-0980
Fax: (905) 764-0982
www.amico.com

HOSPITAL	HOSPITAL
LOCATION	LOCATION
QTY: 1	() TYPE UNITS AS SHOWN / () TYPE UNITS OPPOSITE

A. NURSE CALL MFGR:	MODEL NO.:	DRWG. NO.:
B. MEDICAL GAS MFGR.:	TYPE CONNECTION:	PRICEBOOK-32
C. FINISH:	CEILING HEIGHT:	DRAWN BY:
		CHECKED BY:
		REV. NO.:
		DATE: