




# MAJESTIC SERIES HORIZONTAL HEADWALL SYSTEM

DRAWING # 36

SINGLE TIER with LIGHT TIER and CHASE  
(M/N: HW14-01L-P-M108)



FINISHED FLOOR LINE

VIEW A-A

IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE \_\_\_\_\_

DATE \_\_\_\_\_ PHONE NO. \_\_\_\_\_

TYPE:  
QUANTITY:

SYSTEM DETAILS

SYMBOL	QTY.	DESCRIPTION
HC	2	HORIZONTAL EQUIP. TRACK
SC	1	SERVICE CHASE, 14"
VT	2	VERTICAL EQUIP. TRACK
O	3	GAS. OXYGEN
A	1	GAS. MED AIR
V	4	GAS. VACUUM
E2	6	RECEPTACLE - DUPLEX RED
N2	2	RECEPTACLE - DUPLEX IVORY
SN1	1	SWITCH, S.P.S.T.
SN3	1	SWITCH, 3 POSITION L.V.
NC	1	PROVISION - NURSE CALL 3G
CB	1	PROVISION - CODE BLUE 1G
LV	2	PROVISION - LOW VOLTAGE 1G
MO	1	PROVISION - MONITOR OUT 2G
ETI	1	ELAPSED TIME INDICATOR
R	1	E.T.L. RESET SWITCH
SLM	4	VACUUM SLIDE TRACK MTD.
OL4	1	4' OVERBED LIGHT
LVC*	1	LOW VOLTAGE CONTROLLER
MC	1	MONITOR CHANNEL, TRACK MTD.
S1	1	20"x12" SHELF, TRACK MTD.
B1	1	19"x5"x3" BASKET, TRACK MTD.
AM	1	ANEROID MANOMETER

\*INSTALLED IN OVERBED LIGHT BY AMICO

BED LOCATOR SYSTEM DETAILS

SYMBOL	QTY.	DESCRIPTION
E2	2	RECEPTACLE - DUPLEX RED
N2	2	RECEPTACLE - DUPLEX IVORY
BR	1	RECEPTACLE - BED ONLY
NL	1	NIGHT LIGHT
SC	1	37 PIN CONNECTOR
TL	1	PROVISION - TELEPHONE 1G

**Amico Corporation**  
85 Fulton Hwy, Ontario  
Rehoboth Beach, CANADA  
L4B 2N4, CANADA  
Tel: (905) 877-462-6426 (T)  
Toll-Free: 1-877-462-6426 (T)  
Fax: (905) 764-0960  
www.amico.com

HOSPITAL	HOSPITAL	MODEL NO.:	DRG. NO.
LOCATION	LOCATION	A. NURSE CALL MFGR.:	PRICEBOOK -39
QTY.:	TYPE	B. MEDICAL GAS MFGR.:	DRAWN BY:
1	( ) TYPE UNITS AS SHOWN / ( ) TYPE UNITS OPPOSITE	C. FINISH:	CHECKED BY:
		CEILING HEIGHT:	DATE: