


SELECT SERIES RECESSED WALL CONSOLE

(MIN: E-CON-05)

DRAWING # 50



TYPE:
QUANTITY:

SYSTEM DETAILS	
SYMBOL	DESCRIPTION
O	GAS, OXYGEN
A	GAS, MED AIR
V	GAS, VACUUM
B	BLANK

- GENERAL NOTES:
1. PUNCH OUT CENTER TO CENTER IS 2.50"
 2. THE BACK BOX WILL INCORPORATE 1/2" AND 3/4" KNOCKOUTS ON TOP AND BOTTOM. THE BACK BOX IS PROVIDED WITH ADJUSTABLE FLANGES TO ACCOMMODATE VARIOUS WALL THICKNESS.
 3. FASCIA WILL BE EXTRUDED ANODIZED ALUMINUM.

IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE _____

DATE _____

PHONE NO. _____

A
MICO Corporation
85 Fulton Hwy,
Rehoboth Hill, Ontario
L4B 2N4, CANADA
Tel: (905) 877-462-6426(T)
Toll-Free: 1-877-462-6426(T)
Fax: (905) 764-0960
www.dnisco.com

HOSPITAL _____

LOCATION _____

QTY: () TYPE UNITS AS SHOWN / () TYPE _____ UNITS OPPOSITE

A. NURSE CALL MFGR: _____ MODEL NO.: _____

B. MEDICAL GAS MFGR.: _____ TYPE CONNECTION: _____

C. FINISH: _____ CEILING HEIGHT: _____

DRWG. NO. _____ PRICEBOOK-50 _____

DRAWN BY/CS _____ CHECKED BY/CS _____

REV. NO.: _____ DATE: _____