


# SELECT SERIES RECESSED WALL CONSOLE

DRAWING # 52

(M/N: E-CON-15)



TYPE:  
QUANTITY:

SYSTEM DETAILS	
SYMBOL	DESCRIPTION
E2	RECEPTACLE - DUPLEX RED
NC	PROVISION - NURSE CALL
SLM	VACUUM SLIDE, MOUNTED
SN2	SWITCH, S.P.S.T
O	GAS, OXYGEN
A	GAS, MED AIR
V	GAS, VACUUM
B	BLANK

- GENERAL NOTES:
1. PUNCH OUT CENTER TO CENTER IS 2.50"
  2. THE BACK BOX WILL INCORPORATE 1/2" AND 3/4" KNOCKOUTS ON TOP AND BOTTOM. THE BACK BOX IS PROVIDED WITH ADJUSTABLE FLANGES TO ACCOMMODATE VARIOUS WALL THICKNESS.
  3. FASCIA WILL BE EXTRUDED ANODIZED ALUMINUM.

IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE \_\_\_\_\_ DATE \_\_\_\_\_ PHONE NO. \_\_\_\_\_

**A** *Amico Corporation*  
 85 Fulton Hwy, Ontario  
 Richmond Hill, Canada  
 L4B 2N4, CANADA  
 Toll-Free: 1-877-462-6428 (T)  
 Tel: (905) 764-0800  
 Fax: (905) 764-0882  
 www.amico.com

HOSPITAL \_\_\_\_\_  
 LOCATION \_\_\_\_\_  
 QTY: ( ) TYPE  UNITS AS SHOWN / ( ) TYPE \_\_\_\_\_ UNITS OPPOSITE

A. NURSE CALL MFGR.: \_\_\_\_\_ MODEL NO.: \_\_\_\_\_  
 B. MEDICAL GAS MFGR.: \_\_\_\_\_ TYPE CONNECTION: \_\_\_\_\_  
 C. FINISH: \_\_\_\_\_ CEILING HEIGHT: \_\_\_\_\_  
 DRWG. NO. PRICEBOOK-52  
 DRAWN BY:CS  
 CHECKED BY:CS  
 REV. NO.:  
 DATE: