

FREEDOM SERIES 12" FLOOR PEDESTAL



(M/N: W-PEDESTAL-12)

TYPE:
QUANTITY:

SYSTEM DETAILS	
SYMBOL	DESCRIPTION
PED	FLOOR PEDESTAL, 12"
E2	RECEPTACLE - DUPLEX RED
SLM	VACUUM SLIDE, MOUNTED
O	GAS, OXYGEN
A	GAS, MED AIR
V	GAS, VACUUM

MATERIAL: 304 STAINLESS STEEL (Ga.#16)

FINISH: #4 Satin Finish



IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE _____

DATE _____ PHONE NO. _____

A *Amico Corporation*
 85 Fulton Hwy,
 Richmond Hill, Ontario
 L4B 2N4, CANADA
 Toll-Free: 1-877-462-6426 (T)
 Tel: (905) 764-0900
 Fax: (905) 764-0982
 www.amico.com

HOSPITAL _____
 LOCATION _____
 QTY: (A) TYPE _____ UNITS AS SHOWN / (B) TYPE _____ UNITS OPPOSITE

A. NURSE CALL MFGR: _____ MODEL NO.: _____
 B. MEDICAL GAS MFGR.: _____ TYPE CONNECTION: _____
 C. FINISH: _____ CEILING HEIGHT: _____
 DRWG. NO. _____
 PRICEBOOK-65 _____
 DRAWN BY/CS _____
 CHECKED BY/CS _____
 REV/NO.: _____
 DATE: _____