



# FREEDOM SERIES 18" FLOOR PEDESTAL

(M/N: W-PEDESTAL-18)



TYPE: \_\_\_\_\_  
 QUANTITY: \_\_\_\_\_

SYSTEM DETAILS	
SYMBOL	DESCRIPTION
PED	FLOOR PEDESTAL, 18"
E2	RECEPTACLE - DUPLEX RED
O	GAS, OXYGEN
A	GAS, MED AIR
V	GAS, VACUUM

MATERIAL: 304 STAINLESS STEEL (Ga #16)

FINISH: #4 Satin Finish

IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE \_\_\_\_\_

DATE \_\_\_\_\_

PHONE NO. \_\_\_\_\_

**AMICO** Corporation  
 85 Fulton Way  
 Richmond Hill, Ontario  
 L4B 2N4, CANADA  
 Toll-Free: 1-877-462-6428 (T)  
 Tel: (905) 764-0800  
 Fax: (905) 764-0882  
 www.amico.com

HOSPITAL	A. NURSE CALL MFGR: _____	MODEL NO.:	DRWG. NO.
LOCATION	B. MEDICAL GAS MFGR.: _____	TYPE CONNECTION:	DRAWN BY/CS
QTY: _____	C. FINISH: _____	CEILING HEIGHT: _____	CHECKED BY/CS
( ) TYPE _____ UNITS AS SHOWN / ( ) TYPE _____ UNITS OPPOSITE			REV. NO.:
			DATE: