

- ■ ■ The Amico EMS Digital Oxygen Monitor incorporates the latest microprocessor technology for alarm and surveillance systems. It is designed for use in EMS transportation vehicles and to provide the highest level of user-friendliness and reliability. The digital monitor is constructed to measure gas tank levels and provide an alert when there is a gas deficiency.

- How does it work?

A regulator and transducer are attached to the medical gas cylinder. The transducer sends an electrical signal to the monitor, which displays the pressure and sets off the alarm, if the reading is below the set point.

# EMS



## ■ ■ ■ Key Features

- Microprocessor digital gas pressure display
- Ability to connect the monitor to another unit
- Illuminated LED display readable in any lighting conditions
- Field adjustable alarm set point of 500 psi
- Gas specific regulator with CGA nut and nipple
- Availability in various styles and sizes
- 12 VDC supply (required)