



Alarm Retro-Fit



Did you know that the Ohio Medical Gas Alarm's power supply has potential problems of overheating, which can result in smoke and flame?

Did you know that Amico has a solution to this serious problem? Just ask for the "Alarm Retro-Fit".

Did you know that an Alarm Retro-Fit can be done in less than 30 minutes by an experienced installer?

Did you know that Amico has a 5 year warranty on all our products.

Did you know that Amico has Retro-Fit's for Alarms, Outlets and Nitrogen Control Panels? Visit our web site at www.amico.com for more information.



Orders

Orders can be emailed to : orders@amico.com

or

Faxed to : (905) 764-0862