



Alarm-Manifold Power Supply

Did you know Amico's Power Supply has a built in circuit protection. The circuit protection can be reset by resetting the circuit breaker.



Outlet Extension

Did you know that we have a solution for double drywall? It extends 0.375" more than our regular 1.434" barrel.



Oxy Fill Stations - Survey

Did you know we're thinking of changing the color of our Oxy Fill Station from blue to green to help identify the unit.

Your thoughts would be appreciated. Please email us at marketing@amico.com



Floor Pedestal

Did you know we've redesigned our Floor Pedestal. The top comes off so that the heat can escape, while brazing on the inside.



Orders

Orders can be emailed to : orders@amico.com

or

Faxed to : (905) 764-0862