

specifications

The DISS Dust Caps shall be an Amico ALert-1 Series.

The DISS Dust Caps eliminates dust and dirt build-up on the DISS connections while not in use.

The DISS Dust Cap is made from machined aluminum and black anodized for an aesthetic appearance. Also provided is a 6" long chain that attaches to the outlet to prevent the item from getting lost.



model numbers

O-DIS-CAP-OXY

Used for Oxygen gas only

O-DIS-CAP-GAS

Used for Medical Air, Vacuum, Nitrous Oxide, Carbon Dioxide and Nitrogen gases

O-DIS-CAP-EVA

Used for Evacuation gases, WAGD and AGSS.

Note: This cap is made of plastic material.



project

represented by: