



specifications

The gauges shall be an Amico Alert-1 series.

The gauge shall be with dual scale psi and kPa for pressure, and Hg and kPa for Vacuum.



features

- Size: 2"
- Accuracy: +/- 1.5% of full scale.
ASME B40.1 Grade B
- Range: 0 - 30"Vac [-100 - 0 kPa]
0 - 100 psi [0 - 700 kPa]
0 - 300 psi [0-2000 kPa]
- Temperature: Ambient: -13 °F to 140 °F
Media: Max. 140 °F
- Connection: 1/8" Male NPT Center Back
- Bourdon Tube: Copper alloy
Movement: Copper alloy
Dial: White plastic with stop pin.
Pointer: Black aluminum.
Case: Black steel.
Window: Snap-in acrylic.

model numbers

Vacuum (0-30" Hg) : **V-X-GAUGE-VAC**
 Gas (0-100 psi) : **V-X-GAUGE-100**
 Nitrogen (0-300 psi) : **V-X-GAUGE-300**

project

represented by: