

PLEASE SPECIFY (IF NOT SPECIFIED AMICO WILL PROVIDE DEFAULTS LISTED BELOW):

STUD LENGTH (12'-0"):

STUD SIZE (3/8" TYP.):



STUD GAUGE: #16 TYP.

WALL THICKNESS (5/8" TYP.):

# FREEDOM SERIES RECESSED 24" PANEL WALL

FULL LENGTH GENERAL CARE  
(M/N: FW2448-PR-GEN-P)

DRAWING # 0X




TYPE: \_\_\_\_\_  
 QUANTITY: \_\_\_\_\_

IF MIRRORED UNITS ARE REQUIRED THE CUSTOMER MUST INDICATE THE QUANTITY OF UNITS BELOW AND AMICO WILL PROVIDE A SEPARATE DRAWING SHOWING THE MIRRORED LAY OUT AND QUANTITY

QUANTITY MIRRORED: \_\_\_\_\_

SYSTEM DETAILS	
SYMBOL	DESCRIPTION
PR	1 PANEL WALL
O	X AMICO GAS OUTLET, OXYGEN
A	X AMICO GAS OUTLET, MED AIR
V	X AMICO GAS OUTLET, VACUUM
E2	X RECEPTACLE - DUPLEX RED
N2	X RECEPTACLE - DUPLEX IVORY
TL	X PROVISION - TELEPHONE 1G
	X BLANK PLATE
	2 STUD ASSEMBLY



IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL SIGNATURE \_\_\_\_\_ DATE \_\_\_\_\_ PHONE NO. \_\_\_\_\_

**Amico**  
 85 Fulton Way  
 Richmond Hill, Ontario  
 L4B 2M4, CANADA  
 Toll-Free: 1-877-462-6428(T)  
 Tel: (905) 764-0800  
 Fax: (905) 764-0802  
 www.amico.com

HOSPITAL	HOSPITAL
LOCATION	LOCATION
QTY:	UNITS AS SHOWN

A. NURSE CALL MFRG:	MODEL NO.:	DRWG. NO.
B. MEDICAL GAS MFRG:	TYPE CONNECTION:	MADDY--00
C. FINISH:	CEILING HEIGHT:	DRAWN BY:
		CHECKED BY:
		REV. NO.: 00
		DATE: