

Integrated Patient Lift Pendant Specifications




PLP Features and Benefits:

- Extremely easy to transfer a patient from a wheel chair to bed or stretcher
- Eliminates collision between patient lift ceiling track and ceiling mounted booms
- Continuous charging on track ensures lift is always charged and operable
- 1,000 lbs track weight capacity
- Significantly reduces construction cost of ceiling mounted patient lift tracks

GoLift Specifications:

- **Weight Capacity**
 - up to 400 lbs (181 kg)
 - up to 700 lbs (317 kg)
 - up to 1,000 lbs (453 kg)
- **Lift Motor:** 24 VDC
- **Lift Weight**
 - Under 10 lbs (4.5 kg)
- **Infection Control:**
 - capacitive touch control handle
 - seamless design and finish
 - non-porous rubber emergency strap for easy cleaning
- **Battery Indicator:** visual and audio

Slings



Various sling types available along with the ability to customize to your patients needs. The GoLift slings are ideal for lifts, transfers, limb lifts, toileting and early mobility.

Carry Bar



The carry bar allows for easy sling attachment and quick connection of the hand control.

Capacitive Touch Handle Control



The easy to use capacitive touch handle control is a seamless, clean design with easy connection to the carry bar via magnet.