







GE Healthcare Ordering Guide Reference Chart



Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 <p data-bbox="310 846 457 911">Dash 5000 14 lbs (6.4 kg)*</p>	<p data-bbox="680 659 978 781">N/A (The DASH 5000 has a 0 head compatible plate built into the bottom)</p>	 <p data-bbox="1100 440 1207 464">AHM-050</p>	<p data-bbox="1461 383 1507 407">N/A</p>	 <p data-bbox="1787 659 1850 683">RAC2</p>  <p data-bbox="1692 886 1944 911">Cuffs and Patient Leads</p>
		 <p data-bbox="1104 740 1203 764">SSM-050</p>	<p data-bbox="1461 708 1507 732">N/A</p>	
		 <p data-bbox="1062 1105 1245 1130">TTM-05-450-A00</p>	<p data-bbox="1461 1032 1507 1057">N/A</p>	

*Weight provided to help approximate mounting capacity. For reference only.

GE Healthcare Ordering Guide Reference Chart



Device	Mounting Kit	Mounting Platform	Restrictions	Add-Ons
 <p>RAC2 5 lbs (2.3 kg)*</p>	N/A	 <p>GE-RAC2M-MON-00</p>	Mounts on MRS rail with AHM/SSM	N/A
	 <p>GE-RAC4-SI-00 1.5 lbs (0.68 kg)*</p>	 <p>AHM-050</p>	Mounts below with primary device on AHM/SSM	N/A
 <p>SSM-050</p>				
 <p>Cuffs and Patient Leads</p>	N/A	 <p>100609-MON-BB</p>	Mounts on MRS rail with AHM/SSM	N/A

*Weight provided to help approximate mounting capacity. For reference only.