

## Clarity Series

# Clarity 40

Task & Reading Lights



### Features:

**Fitted With:** 1 x incandescent lamp 40 W

**Work Equipment:** Integrated in spare lamp

**Connected Load:** 120 V; 50/60 Hz

**Mains Lead: Approx.** 157" (4 m) cordset

**Luminaire Body:** Polycarbonate (PC)

**Weight (NET):** Approx. 3.1 lbs (1.4 kg)

**Power Consumption:** Approx. 41 W

**Usage:** Rocker switch (I/O)

**Technology:** Single-pole primary

**Class Of Protection:** II

**Norms:** EN 60598-1, EN 60598-2-25

**Special Features:** Separate light strip for orientation

**Fastening:** Rail mount included for rail mounted light

The cost-effective Clarity 40, featuring a 40-Watt incandescent bulb, is cool to the touch even after hours of operation. Use the Clarity 40 for patient room illumination. The light is ideal for reading and provides a special illuminated handle for night time security. The Clarity 40 is cool to the touch after hours of operation. It also incorporates a night light and features an incandescent 40-watt lamp. A 13-watt compact fluorescent lamp can be used with diffuser models.

# Dimensions



## Model Examples

<b>L-CL40-AB-RH-A</b>	Clarity 40, Articulating Arm, Basic, Rail Mount, Horizontal, Amico
<b>L-CL40-AB-RV-F</b>	Clarity 40, Articulating Arm, Basic, Rail Mount, Vertical, Fairfield
<b>L-CL40-EB-RH-A</b>	Clarity 40, Articulating Arm with Extension, Basic, Rail Mount, Horizontal, Amico
<b>L-CL40-EB-RV-F</b>	Clarity 40, Articulating Arm with Extension, Basic, Rail Mount, Vertical, Fairfield

**NOTE:** North American: 120 V; International: 230 V

## Ordering Guide

