

Hydraulic Patient Transfer Stretcher (S-H-300)



The Amico Hydraulic Patient Transfer Stretcher (S-H-300) provides caregivers the control they need in the emergency environment. The stretchers can be adjusted to various positions – raising, lowering, trendelenburg, and reverse trendelenburg – by actuating the hydraulic pistons. The back section is adjustable through the use of gas springs and the knee gatch can be raised and lowered with an optional manual crank.

The stretcher is made of carbon steel and painted by an electrostatic spray (powder-coat). The base cover is made from injectionmolded ABS plastic and the locking siderails are made from a combination of steel, aluminum extrusion and plastic parts for increased durability and to facilitate easy cleaning. The back section includes an X-Ray-compatible board.

The stretcher comes complete with an IV pole, bumper wheels and actuator dust covers to protect the hydraulic actuators.

Key Features:

- Centrally controlled brake system
- Collapsible siderails
- Dual side foot pedals
- Retractable fifth wheel steering system
- X-Ray compatible backrest
- Shock absorbing, non-marking bumpers
- Foley bag holder
- 4" (10 cm) mattress included
- Ergonomic push handles
- Storage area in base for oxygen E-cylinder tank
- Storage tray for patient belongings
- Permanent IV Pole

Model Number & Approval

Model Number	S-H-300
--------------	---------

Overall Dimensions

Overall Length	84.5" (214 cm)
Overall Width	33" (84 cm)

Litter Positioning

Trendelenburg/Reverse Trendelenburg	16°
Maximum Fowler Section Inclination	80° ± 5
Maximum Knee Section Inclination	35° ± 5

Height Range

Lowest	23" (58 cm)
Highest	35" (89 cm)

Surface Dimensions

Surface Width	25" (64 cm)
Surface Length	75" (191 cm)
Surface Thickness	4" (10 cm)

Weight Capacity & Caster Dimension

Maximum Weight Limit	550 lbs (250 kg)
Caster Diameter	8" (20 cm)

Optional Features & Accessories:

- Utility Tray
- Four Sided Brake
- Overbed tables
- Removable IV Pole

Amico Beds Corporation reserves the right to change product specifications without notice.

Warranty: Three years parts

Amico Beds Corporation 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada | 600 Prime Place, Hauppauge, NY 11788, USA
Toll Free Tel: 1.877.462.6426 | Toll Free Fax: 1.866.440.4986 | Tel: 905.764.0800 | Fax: 905.764.0862