

Teaching & Training Solutions

Sapphire/Regal Surface Mount Headwalls

Product Benefits:

Amico's designers have developed a revolutionary headwall that is self-contained and alleviates the headaches associated with a nursing simulation lab installation.

The traditional method of outfitting a simulation lab with medical gas services can cost thousands of dollars for pumps, compressors and a network of plumbing. This Plug & Play training headwall simulates Oxygen, Air and even Vacuum just by connecting it to a portable compressor.


The compressor can be located anywhere you wish, you can even provide your own in-wall source. It is the quickest and most inexpensive method for simulating a patient room's medical devices.




Product Applications:

Educational Facilities
Teaching/Training
Universities
Simulation Centers
Vocational Schools
Nursing Labs
Mock Up Rooms
Institutions

Standard Colors:



Antique White



Fusion Maple



Wild Cherry

Features:

- Plug & Play headwall
- Typical widths: 24" (60.96 cm), 36" (91.44 cm), 48" (121.92 cm), 60" (152.4 cm)
- One full length, scratch-resistant accessory rail
- Three standard colors of hospital grade laminate to choose from
- Super silent, HP air compressor to simulate medical gases
- Rear connection option available for electrical and/or compressor and vacuum
- Mobile floor stand option
- Optional light for units 48" or wider


Distributed By:




Teaching & Training Solutions

Sapphire/Regal Surface Mount Headwalls

Sapphire Teaching/Training Headwall



Regal Teaching/Training Flatwall



Sapphire with light Teaching/Training Headwall



Item	Description	Item	Description
1	Duplex Receptacle	5	Nurse Call Cover
2	Vacuum Switch	6	Compressed Air Outlet
3	Vacuum Outlet	7	Medical Air Outlet
4	Oxygen Outlet	8	Light (Optional)

*Not shown